



655 Research Parkway
Suite 525
Oklahoma City, OK 73104
Phone: 405-271-1544
Fax: 405-271-1545
E-mail: info@hyalose.com
www.hyalose.com

Specification Sheet

Select-HA™ LoLadder

Catalogue No. HYA-LOLAD-20

Product Description

Select-HA™ is a hyaluronic acid (HA) preparation of uniform and narrow size distribution prepared by *in vitro* synthesis using recombinant *Pasteurella multocida* hyaluronan synthase¹. The LoLadder contains five Select-HA™ molecular mass markers in the range of ~25 kDa to 500kDa.

Quantity

The tube contains 25µg (total) of five different Select-HA™ preparations supplied as a mixture of lyophilized sodium salts. The materials are sufficient for at least 20 lanes.

Recommended Reconstitution Procedure

Centrifuge the tube for a few seconds to collect the Select-HA™ solids in the bottom of the tube. Carefully open and add 100 µl of sterile water directly to the bottom of the tube. Allow two hours at 4°C for sample rehydration and then mix well before use. Electrophoresis of 5µl of reconstituted sample (Fig. 1) using a standard gel loading buffer on an agarose gel (0.6 – 1%) results in clearly defined bands when stained with 0.005% Stains-all (in 50% ethanol)².

Storage

Store at -20°C or below. To avoid frequent freeze-thaw, aliquoting is recommended. Avoid contamination with microbes or HA-degrading enzymes.

References

1. DeAngelis, P.L. *et al.* (1998) *J. Biol. Chem.*, **273**, 8454–8458.
2. Lee, H.G. *et al.* (1994) *Analytical Biochem.* **219**, 278-287.

For research use as molecular mass markers only

U.S. patent 6,444,447 and patents pending